

THE NATIONAL INSTITUTE OF ENGINEERING MYSURU-570006

An Autonomous Institute Affiliated to VTU, Belagavi



Dissertation On

NEURO PREDICT

Submitted in partial fulfillment for the award of the degree of Bachelor of Engineering in
Electronics and Communication Engineering

Submitted by

ANNANGI SAIKIRAN BABU USN:4NI19EC124


K P YASHIKA USN:4NI19EC127

K S NISHAAN KUSHALAPPA USN:4NI19EC044

ANIRUDDH ADIGA USN:4NI19EC014

Under the guidance of
DR. REMYA JAYACHANDRAN

Assistant Professor
Department of Electronics and Communication



DEPARTMENT OF ELECTRONICS AND COMMUNICATION ENGINEERING

2022-23

THE NATIONAL INSTITUTE OF ENGINEERING MYSURU-570006

An Autonomous Institute Affiliated to VTU, Belagavi

CERTIFICATE

This is to certify that the project entitled “NeuroPredict” has been successfully completed by **ANNANGI SAIKIRAN BABU**(4NI19EC124), **K P YASHIKA**(4NI19EC127), **K S NISHAAN KUSHALAPPA**(4NI19EC044), **ANIRUDDH ADIGA**(4NI19EC014) of 8th semester B.E. who carried out the project work under the guidance of **DR.REMYA JAYACHANDRAN** and **DR. NARASIMHA KAULGUD** in the Partial fulfillment for the award of the degree of Bachelor Engineering in Electronics and Communication Engineering of Visvesvaraya Technological University, Belagavi during the year 2022-23. It is certified that all corrections/Suggestions indicated during the internal assessment have been incorporated into the report. The Project has been approved in partial fulfillment for the award of the said degree as per academic regulations of the National Institute of Engineering (An Autonomous Institution under Visvesvaraya Technological University, Belagavi).

Dr. Remya Jayachandran

Project Guide

Department of ECE

Dr. S Parameshwara

Head of the Department

Department of ECE

Dr. Rohini Nagapadma

Principal

NIE, Mysuru

External Viva-

Sl. No.	Name of the Examiners	Signature
1		
2		
3		

DECLARATION

We, **ANNANGI SAIKIRAN BABU**(4NI19EC124), **K P YASHIKA**(4NI19EC127), **NISHAAN KUSHALAPPA**(4NI19EC044), **ANIRUDDH ADIGA**(4NI19EC014) students of 8th semester of Bachelor of Engineering, Department of Electronics and Communication, The National Institute of Engineering, Mysuru hereby declare that the project work entitled "**"NUERO-PREDICT"**" has been carried out by us under the guidance of Dr. Remya Jayachandran (Assistant Professor, Department of ECE). This project work is submitted to The National Institute of Engineering, Mysuru, (An Autonomous institute under VTU, Belagavi) in partial fulfillment of the course requirements for the award of degree in Bachelor of Engineering during the academic year 2022-23. This written submission represents a record of original work and we have adequately cited and referenced the original sources.

Furthermore, the matter embodied in this thesis has not been submitted to any other University or Institution for the award of any degree.

Place: Mysuru

Date:

(Signature of the Students)

ACKNOWLEDGEMENT

It gives immense pleasure in bringing out this synopsis of the project entitled “Neuro Predict” major project report as part of the curriculum of Electronics and Communication Engineering. I wish to thank all the people who gave me unending support.

Firstly, we would like to thank our teacher and guide **Dr. Remya Jayachandran** and **Dr Narasimha Kaulgud**, Department of Electronics and Communication, NIE, Mysuru, who gave us her valuable suggestions and ideas when we were in need of them. They encouraged us to work on this project. We wish to express our sincere gratitude to our project guide for constant guidance and encouragement in every stage of preparation. Their careful monitoring of the complete process helped in improving in presentation of the ideas. We take this opportunity to sincerely acknowledge her for performing all evaluation studies.

We would like to express our gratitude towards **Dr. Rohini Nagapadma**, The Principal, NIE for extending all necessary facilities.

We would like to express my gratitude towards **Dr. S Parameshwara**, Head of the Department, ECE, NIE for providing all necessary facilities.

Finally, we thank all my professors, non-teaching staff, and my parents, for their constant support and guidance. We are immensely grateful to all involved in this project as without their inspiration and valuable suggestion it would not have been possible to develop the project within the prescribed time.

With utmost gratitude,

Annangi Saikiran Babu

K P Yashika

K S Nishaan Kushalappa

Aniruddh Adiga

ABSTRACT

Alzheimer's disease is a global health concern requiring early prediction and diagnosis for better patient outcomes. The proposed NeuroPredict system utilizes a Convolutional Neural Network (CNN) model to accurately predict Alzheimer's disease based on MRI images. Leveraging CNN's image recognition capabilities, NeuroPredict offers a robust solution for efficient and accurate prediction. Exploratory Data Analysis (EDA) is emphasized to uncover patterns and relationships within the Alzheimer patient's dataset, providing insights for research and interventions. The project also introduces a hardware kit for proactive prevention and personalized support, along with a user-friendly web application for seamless Alzheimer's prediction. NeuroPredict aims to impact Alzheimer's care and advance neurodegenerative disease research, contributing to improved management and outcomes.

TABLE OF CONTENTS

1	INTRODUCTION	4
2	LITERATURE SURVEY	6
2.1	Early prediction of Alzheimer's disease using convolutional neural network: a review.	6
2.2	Alzheimer's Disease Neuroimaging Initiative. Convolutional Neural Network Based MRI Image Analysis for the Alzheimer's Disease Prediction From Mild Cognitive Impairment.	6
2.3	Alzheimer's Disease Prediction Using Convolutional Neural Network Models Leveraging Pre-existing Architecture and Transfer Learning	7
2.4	Convolutional Neural Networks for Alzheimer's Disease Detection on MRI Images	8
3	Difference in MRI scans	9
4	CNN-Convolutional neural network	10
4.1	Dataset description	10
4.2	Data Preprocessing	12
4.3	Model Building	13
4.4	Model Parameters	16
4.5	Training and Evaluation	17
4.6	CNN Model accuracy and results	19
5	Exploratory Data-analysis on the Alzheimer's dataset	21
5.1	Problem Statement	21
5.2	Need of exploratory data analysis on Alzheimer's dataset	21

5.3	Data Understanding	22
5.4	Dataset details	22
5.5	Stages of performing Exploratory data analysis	23
5.6	Distribution of groups	24
5.7	Distribution of Clinical Dementia Rating	24
5.8	Mini-Mental State Examination	26
5.9	Age-Group	27
5.10	Bivariate Analysis	27
5.11	Clinical Dementia rating	28
5.12	Years of education	28
5.13	Social economic status	29
5.14	Number of visits	31
5.15	Summary of Observations:	32
6	Website	34
6.1	Website Overview	34
6.2	User Interface Design	36
6.3	Functionality and Features	36
6.4	Website architecture and technologies	37
6.5	Conclusion	38
7	Wearable Devices - Medical assisting kit	41
7.1	Problem Statement	42
7.2	Objective	42
7.3	Components required	43

7.4	Methodology	43
7.5	MemorEYEs: Wearable Device for Alzheimer's Patients	44
7.5.1	Smart Spectacles	44
7.5.2	Wrist Band	46
7.6	Technological Considerations	46
7.7	Advantages of MemorEYEs	48
7.8	Future Scope- Medical assisting kit	49
7.9	Integration with Data analysis	50
7.10	Conclusion- Medical assisting kits	51
8	Future Scope-NeuroPredict	52
9	Conclusions - NeuroPredict	53
	REFERENCES	55

LIST OF FIGURES

3.1	Figure A - MRI scan of a patient with Alzheimer's disease, Figure B - MRI scan of a healthy individual	10
4.1	MRI images of various patients used in training the model	12
4.2	CNN Architecture	15
4.3	Training the CNN Model	18
5.1	Distribution of group variable	25
5.2	Clinical Dementia Rating (CDR) Score interpretations	26
5.3	Stacked bar plot of Mini-Mental State Examination	27
5.4	Stacked bar plot for Age-Group analysis	28
5.5	CDR and Age analysis	29
5.6	CDR plot with EDUC	30
5.7	Stacked bar plot of Social and Economic status	31
5.8	Plot of Social and Economic status	32
5.9	Number of visits	33
6.1	Uploading MRI image.	39
6.2	Output shown for the MRI image uploaded.	41
7.1	The designed Smart Spectacles prototype	45
7.2	Block diagram of Wrist Band	47
7.3	The designed Wrist Band prototype	48

1 INTRODUCTION

Alzheimer's disease is a significant global health concern, and early prediction and diagnosis play a crucial role in improving patient outcomes. To address this challenge, our project, "NeuroPredict", utilizes a Convolutional Neural Network (CNN) model to predict Alzheimer's disease based on MRI images. CNNs have demonstrated exceptional performance in image recognition tasks, making them well-suited for analyzing medical images like MRI scans. By leveraging the power of CNN, NeuroPredict offers a robust solution for accurate and efficient prediction of Alzheimer's disease.

In addition to the CNN model, our project also emphasizes the importance of Exploratory Data Analysis (EDA) in understanding the underlying patterns and relationships within the Alzheimer's dataset. Through EDA, we delve into demographics, clinical information, and derived anatomical volumes to gain insights into the factors contributing to the development and progression of the disease. By identifying correlations and uncovering meaningful information, EDA provides a foundation for further research, decision-making, and the development of targeted interventions.

The significance of this project lies in the far-reaching impact of Alzheimer's disease. Timely prediction and diagnosis of the disease are essential for enabling early interventions, treatment planning, and support for affected individuals and their families. By harnessing the potential of NeuroPredict and its accurate prediction capabilities, we aim to improve early detection and diagnosis, ultimately leading to better management and outcomes for Alzheimer's patients.

Furthermore, our project goes beyond prediction and diagnosis. We have also developed a hardware kit that utilizes the insights obtained from Exploratory Data Analysis (EDA). This kit serves as a proactive approach to prevent Alzheimer's in vulnerable groups and assist individuals already diagnosed with the disease. By analyzing the relationships between various parameters such as gender, years of education, socioeconomic status, and age, our Exploratory Data Analysis (EDA) has enabled the development of a comprehensive solution. The hardware kit targets preventive measures and provides tailored support to individuals based on their specific needs, contributing to a more proactive and personalized approach to Alzheimer's care.

To ensure accessibility and convenience, we have also created a user-friendly web appli-

cation as part of our project NeuroPredict. This web application allows users to directly upload their MRI images, making the Alzheimer's prediction process seamless and efficient. By leveraging the power of advanced technology, data analysis, and user-centric design, we have made Alzheimer's prediction accessible to individuals, empowering them to make informed decisions about their health and seek appropriate medical care.

In conclusion, our project combines the capabilities of CNN models, the insights obtained from EDA, and the development of innovative hardware and web applications. By focusing on early prediction and diagnosis, proactive prevention, and personalized care, our project aims to make a significant impact on the lives of individuals affected by Alzheimer's disease. Moreover, it contributes to the broader research and understanding of neurodegenerative diseases, driving advancements in the field and improving the overall healthcare landscape.

2 LITERATURE SURVEY

2.1 Early prediction of Alzheimer's disease using convolutional neural network: a review.

[1] Patil, V., Madgi, M. Kiran, A. Early prediction of Alzheimer's disease using convolutional neural network: a review. Egypt J Neurol Psychiatry Neurosurg 58, 130 (2022). <https://doi.org/10.1186/s41983-022-00571-w>

This paper provides a comprehensive review of convolutional neural networks (CNNs) for early Alzheimer's disease (AD) prediction. The authors discuss the limitations of traditional diagnostic methods and highlight CNNs as a promising approach in AD diagnosis. The review covers various neuroimaging modalities and CNN architectures used in previous studies. It summarizes performance metrics such as accuracy, sensitivity, specificity, and AUC, demonstrating the superior performance of CNN models over traditional methods. The paper also addresses challenges in this field, including limited datasets and model interpretability, and suggests future research directions. It emphasizes the importance of larger datasets and the potential of hybrid models integrating clinical, genetic, and neuroimaging data.

2.2 Alzheimer's Disease Neuroimaging Initiative. Convolutional Neural Network Based MRI Image Analysis for the Alzheimer's Disease Prediction From Mild Cognitive Impairment.

[2] Lin W, Tong T, Gao Q, Guo D, Du X, Yang Y, Guo G, Xiao M, Du M, Qu X; Alzheimer's Disease Neuroimaging Initiative. Convolutional Neural Networks-Based MRI Image Analysis for the Alzheimer's Disease Prediction From Mild Cognitive Impairment. Front Neurosci. 2018 Nov 5;12:777. doi: 10.3389/fnins.2018.00777. PMID: 30455622; PMCID: PMC6231297.

This paper explores the use of Convolutional Neural Networks (CNNs) to predict Alzheimer's disease (AD) progression from mild cognitive impairment (MCI) using MRI images. The study

utilized data from the Alzheimer's Disease Neuroimaging Initiative (ADNI) and employed a CNN-based framework for feature extraction and classification. Results demonstrated the effectiveness of the CNN approach in predicting AD, outperforming traditional machine learning algorithms. The study's strengths include the large dataset and the automated nature of the CNN-based framework. However, limitations include the exclusive reliance on MRI images and the need for external validation. Overall, the paper contributes to the application of deep learning in neuroimaging for early AD detection, showing promising results that require further validation.

2.3 Alzheimer's Disease Prediction Using Convolutional Neural Network Models Leveraging Pre-existing Architecture and Transfer Learning

[3] M. T. Abed, U. Fatema, S. A. Nabil, M. A. Alam and M. T. Reza, "Alzheimer's Disease Prediction Using Convolutional Neural Network Models Leveraging Pre-existing Architecture and Transfer Learning," 2020 Joint 9th International Conference on Informatics, Electronics Vision (ICIEV) and 2020 4th International Conference on Imaging, Vision Pattern Recognition (icIVPR), Kitakyushu, Japan, 2020, pp. 1-6, doi: 10.1109/ICIEVicIVPR48672.2020.9306649.

This paper explores the application of convolutional neural network (CNN) models for Alzheimer's disease prediction. Their research focuses on leveraging pre-existing architectures and transfer learning to improve prediction accuracy. The study demonstrates the effectiveness of CNN models in analyzing medical imaging data for Alzheimer's disease prediction. By extracting meaningful features from brain images and training a predictive model, the authors achieve promising results in terms of accuracy and performance. The experiments validate the potential of CNN models, in combination with transfer learning, to enhance prediction tasks. The research contributes to the field of medical image analysis and provides valuable insights for developing robust predictive models for Alzheimer's disease. Overall, the paper highlights the significance of CNN models and transfer learning techniques in Alzheimer's disease prediction, offering valuable insights for researchers and practitioners in the field.

2.4 Convolutional Neural Networks for Alzheimer's Disease Detection on MRI Images

[4] Ebrahimi A, Luo S; Alzheimer's Disease Neuroimaging Initiative. Convolutional neural networks for Alzheimer's disease detection on MRI images. *J Med Imaging* (Bellingham). 2021 Mar;8(2):024503. doi: 10.1117/1.JMI.8.2.024503. Epub 2021 Apr 29. PMID: 33937437; PMCID: PMC8083897.

This paper investigates the application of convolutional neural networks (CNNs) for detecting Alzheimer's disease using MRI images. The study focuses on developing a CNN-based approach to accurately classify Alzheimer's disease and healthy control subjects. The authors utilize a large dataset from the Alzheimer's Disease Neuroimaging Initiative (ADNI) and demonstrate the effectiveness of CNNs in distinguishing between Alzheimer's disease and normal brain images. Their experimental results show promising performance, with high accuracy achieved in the classification task. The paper highlights the potential of CNNs as a powerful tool for early detection and diagnosis of Alzheimer's disease, which can significantly aid in timely interventions and treatment. The research contributes to the field of medical imaging and presents a valuable approach for improving Alzheimer's disease detection using MRI images. Overall, this paper provides important insights into the application of CNNs for Alzheimer's disease detection, paving the way for further advancements in this area of research.

3 Difference in MRI scans

MRI (Magnetic Resonance Imaging) scans can reveal distinct differences between the brain of a person with Alzheimer's disease (AD) and a healthy individual. These differences are critical in the accurate diagnosis and understanding of the pathological changes associated with AD. Here are some key differences typically observed in MRI scans:

Hippocampal Atrophy: One of the hallmark signs of AD is the shrinkage or atrophy of the hippocampus, a region in the brain involved in memory formation. MRI scans of individuals with AD often show significant volume loss or reduced size of the hippocampus compared to healthy individuals.

Cortical Thinning: AD affects the cerebral cortex, the outer layer of the brain responsible for various cognitive functions. MRI scans can detect cortical thinning, particularly in regions associated with memory, language, and executive functions. This thinning is indicative of the degenerative changes occurring in the brain due to AD.

Ventricular Enlargement: As brain tissue degenerates in AD, the ventricles (fluid-filled spaces within the brain) tend to enlarge. This expansion is visible in MRI scans and is linked to the loss of surrounding brain tissue.

White Matter Abnormalities: White matter, composed of nerve fibers that connect different regions of the brain, can be affected by AD. MRI scans may reveal white matter hyperintensities or abnormal signal intensities, indicating damage or disruption in the integrity of these fiber pathways.

Amyloid Plaques and Tau Tangles: While MRI scans cannot directly detect amyloid plaques and tau tangles, the characteristic protein abnormalities in AD, advanced imaging techniques like amyloid PET scans or tau PET scans can be used in conjunction with MRI to provide a more comprehensive assessment of AD pathology.

It is important to note that these MRI differences are observed at the group level, and individual variations exist. Additionally, these MRI findings are not specific to AD and can be seen in other neurodegenerative conditions as well. Therefore, a comprehensive diagnosis of AD

typically involves combining MRI findings with clinical evaluations, cognitive assessments, and other diagnostic tests.

Overall, MRI plays a vital role in identifying structural changes in the brains of individuals with AD, aiding in early detection, monitoring disease progression, and informing treatment strategies. These MRI differences provide valuable insights into the underlying neurodegenerative processes associated with AD, contributing to a better understanding of the disease and the development of potential therapeutic interventions.




Figure 3.1: Figure A - MRI scan of a patient with Alzheimer's disease,
Figure B - MRI scan of a healthy individual

4 CNN-Convolutional neural network

4.1 Dataset description

Dementia, a progressive brain disorder characterized by memory loss and cognitive decline, is a significant public health concern. It is the sixth leading cause of death in the United States, and millions of Americans are affected by this condition. While there is currently no cure for dementia, early diagnosis and intervention can help slow its progression and improve the quality of life for individuals living with the disease. Magnetic resonance imaging (MRI) is a non-invasive medical

test that utilizes a strong magnetic field and radio waves to generate detailed images of the brain. It plays a crucial role in the diagnosis and monitoring of dementia. MRI scans can provide valuable insights into the structural and functional changes that occur in the brain, aiding in the identification and characterization of dementia subtypes. To leverage the power of machine learning, you have obtained access to a dataset from the Alzheimer's Disease Neuroimaging Initiative (ADNI). This dataset consists of MRI images of individuals both with and without dementia, and each image is labeled with the corresponding diagnosis. By utilizing this data, you can train a deep learning model to classify MRI images and potentially develop an automated system for dementia diagnosis.

The dataset is divided into two sets: a training set and a testing set. The training set contains a significant number of images, including 717 mild demented, 52 moderate demented, 2560 non-demented, and 1792 very mild demented images. The testing set, on the other hand, consists of 179 mild demented, 12 moderate demented, 640 non-demented, and 448 very mild demented images. This division allows you to train your model on a large and diverse set of images and evaluate its performance on unseen data. Deep learning, a subset of machine learning, employs artificial neural networks to learn complex patterns and representations directly from data. These networks, often referred to as deep neural networks, are capable of automatically extracting relevant features from the input images. Convolutional neural networks (CNNs) are particularly effective for image classification tasks, making them a suitable choice for analyzing MRI images in the context of dementia. Before feeding the MRI images into the deep learning model, it is necessary to preprocess the data. Preprocessing steps may involve resizing the images to a consistent resolution, normalizing pixel values to a common range, and applying data augmentation techniques to enhance the diversity and generalization of the dataset. These steps help optimize the model's learning process and improve its ability to generalize to new data. Once the dataset is preprocessed, you can select an appropriate deep learning model architecture. Popular choices for image classification tasks include CNN architectures like AlexNet, VGGNet, ResNet, Inception, or DenseNet. Depending on the complexity and size of your dataset, you may opt for a simpler or more advanced architecture. Experimenting with different architectures and hyperparameters can help you find the optimal model configuration for your specific task.

With the dataset split and the model architecture chosen, you can proceed to train the deep learning model. Training involves presenting the labeled MRI images from the training set to the

model iteratively, adjusting the model's internal parameters through a process called backpropagation. The goal is for the model to learn and discriminate the features that differentiate between dementia and non-dementia cases, enabling accurate classification. After training the model, it's crucial to evaluate its performance on the testing dataset. This evaluation provides insights into how well the model generalizes to new, unseen data. Metrics such as accuracy, precision, recall, and F1 score can be used to assess the model's performance and compare it against other approaches. Further analysis, such as generating confusion matrices, can help identify potential areas of improvement and guide future research.




Figure 4.1: MRI images of various patients used in training the model

4.2 Data Preprocessing

Data preprocessing plays a crucial role in the data analysis pipeline, as it involves a series of steps to clean, transform, and format raw data before it can be used for analysis. This process is essential to ensure the accuracy, consistency, and completeness of the data, as well as to mitigate any potential biases or errors that may arise from the initial data collection. Data cleaning is the first and fundamental step in data preprocessing. It involves identifying and rectifying errors, inconsistencies, and outliers within the dataset. This process can be time-consuming and challenging, requiring careful examination and verification of each data point. Common techniques used in data cleaning include identifying and removing outliers, imputing missing values, and correcting errors. Outliers are data points that significantly deviate from the rest of the dataset and can have a disproportionate impact

on analysis results. By identifying and removing these outliers, the overall integrity and reliability of the data are enhanced. Imputing missing values is another important aspect of data cleaning, as missing data can hinder the analysis process. Various methods, such as mean imputation, regression imputation, or sophisticated machine learning algorithms, can be employed to estimate and fill in missing values based on the available information. Additionally, correcting errors in the data is crucial to ensure its accuracy. Errors can arise due to various factors, such as human input mistakes or technical issues during data collection. By carefully examining the data for inconsistencies and errors and applying appropriate corrective measures, the reliability of subsequent analyses can be improved. After data cleaning, the dataset may require transformation to make it suitable for analysis. This involves converting data types, scaling values, or normalizing data to ensure comparability and compatibility across different features and variables. Data type conversion is necessary to represent the data in a format that aligns with the analysis requirements. For example, converting text-based data into numerical or categorical formats allows for mathematical operations or statistical analysis. Scaling data involves adjusting the range of values within a feature to a common scale, which can be helpful when different features have significantly different value ranges. Normalizing data adjusts the distribution of values to have a mean of 0 and a standard deviation of 1, which can be beneficial when working with algorithms that assume normalized data or when comparing variables with different units. Finally, formatting the preprocessed data in a suitable structure is essential for effective analysis. This may involve converting the data into a tabular format, spreadsheet, or database, depending on the tools and techniques to be employed. Ensuring the data is well-organized and easily accessible facilitates subsequent analyses and allows for seamless integration with various statistical and machine learning models. Data preprocessing is a vital step in the data analysis pipeline, as it sets the foundation for accurate and reliable insights. Through careful data cleaning, transformation, and formatting, potential biases, errors, and inconsistencies can be mitigated, enabling researchers and analysts to derive meaningful conclusions and make informed decisions. In conclusion, data preprocessing encompasses a series of tasks that ensure the data is clean, transformed, and formatted appropriately for analysis purposes. By paying attention to each step of the preprocessing pipeline, researchers can enhance the quality of their analyses and uncover valuable insights from the data.

4.3 Model Building

The "cnn-model" is a well-structured and powerful convolutional neural network (CNN) architecture specifically designed for image classification tasks. It follows a sequential design pattern and comprises a total of 10 layers, including convolutional layers, pooling layers, and fully connected layers. Convolutional layers are the backbone of the model, responsible for extracting relevant features from the input image. These layers employ a series of filters or kernels that convolve across the image, enabling the network to detect patterns, edges, and textures at various scales and orientations. Pooling layers play a vital role in reducing the spatial dimensions of the feature maps generated by the convolutional layers. By downsampling the feature maps, pooling layers help to capture the most salient information while reducing computational complexity and controlling overfitting. Common pooling techniques include max pooling and average pooling. The fully connected layers follow the convolutional and pooling layers and serve as the classification head of the model. These layers take the high-level feature representations learned by the preceding layers and map them to specific output classes. By leveraging these layers, the model can understand and distinguish between different objects or concepts within the input image. The "cnn-model" architecture also incorporates appropriate activation functions to introduce non-linearities into the network. For the convolutional layers, the "act" activation function is utilized. This non-linear activation function enhances the model's ability to capture intricate and complex patterns present in the image data. To ensure accurate classification probabilities, the output layer of the model employs the "softmax" activation function. This specialized activation function produces a probability distribution across the different output classes, enabling the model to assign a likelihood score to each class. In the "cnn-model" architecture, the "conv-block" function is employed to define the convolutional layers. This function takes the number of filters, the kernel size, and the padding as inputs, allowing flexibility in specifying the desired characteristics of each convolutional layer. The "act" activation function is applied within this function to introduce non-linearity and increase the model's capacity to learn intricate features.

Similarly, the "dense-block" function is used to define the fully connected layers. This function takes parameters such as the number of neurons, the dropout rate, and the activation function. Dropout is an important regularization technique that randomly sets a fraction of input units




Figure 4.2: CNN Architecture

to zero during training, reducing the risk of overfitting and enhancing the model's generalization ability.

To train the "cnn-model," the widely-used "fit" method from the "keras.models.Model" class is employed. This method takes the training data, training labels, the number of training epochs, and the batch size as inputs. By iteratively adjusting the model's internal parameters based on the training data, the model learns to recognize and classify the patterns associated with the given image classes. Evaluation of the trained model is performed using the "evaluate" method of the "keras.models.Model" class. This method takes the test data, test labels, and the batch size as inputs and returns evaluation metrics such as loss and accuracy. It allows for the assessment of the model's performance on unseen data, providing insights into its generalization capabilities. Once trained, the "cnn-model" can be utilized for image classification by invoking the "predict" method of the "keras.models.Model" class. This method takes the input image as input and produces the probability distribution of each class as the output. This enables the model to assign a likelihood score to each class, indicating the predicted probability of the input image belonging to a particular category. In summary, the "cnn-model" is a sophisticated and effective convolutional

neural network (CNN) architecture specifically designed for image classification tasks. With its well-structured layers, including convolutional, pooling, and fully connected layers, the model excels at extracting relevant features from input images and accurately classifying them. By utilizing the "act" activation function for convolutional layers and the "softmax" activation function for the output layer, the model introduces non-linearities and ensures reliable classification probabilities for each output class. The "cnn-model" incorporates the "conv-block" and "dense-block" functions to define the convolutional and fully connected layers, respectively. These functions allow for flexibility in specifying the number of filters, kernel size, padding, number of neurons, dropout rate, and activation function, empowering the model with the ability to learn complex patterns and make informed predictions. Training the "cnn-model" is achieved through the "fit" method, where the model learns from the training data by iteratively adjusting its internal parameters. Evaluation of the trained model is performed using the "evaluate" method, providing valuable metrics such as loss and accuracy on unseen test data. Ultimately, the "cnn-model" is a valuable tool for image classification tasks, capable of accurately predicting the probability of an input image belonging to different classes. Its comprehensive architecture, activation functions, and training capabilities make it an excellent choice for researchers and practitioners in the field of image analysis and classification.

4.4 Model Parameters

The provided model architecture showcases a well-designed framework for image classification. Let's delve deeper into its components and their functionalities. The model commences with an input layer configured to handle grayscale images with dimensions [None, 48, 48]. This specification implies that the model expects input images of size 48x48 pixels. Following the input layer, the model incorporates a Conv2D layer, which employs 32 filters with a kernel size of 3x3. The activation function used in this layer is the popular 'relu' (Rectified Linear Unit), known for its ability to introduce non-linearity and capture complex features within the data. Continuing the architecture, a MaxPooling2D layer is introduced, employing a pool size of 2x2. Max pooling reduces the spatial dimensions of the input, allowing for the extraction of dominant features while retaining important information. Next, another Conv2D layer is added to the model. This time, it

employs 64 filters with a kernel size of 3x3, again utilizing the 'relu' activation function to capture higher-level features. Subsequently, a second MaxPooling2D layer is incorporated, following the same principles as the previous pooling layer to further downsample the feature maps. To prepare the data for classification, a Flatten layer is included, which transforms the output from the previous pooling layer into a one-dimensional vector. Continuing with the architecture, the model introduces a Dense layer with 128 neurons. Dense layers are fully connected layers, where each neuron is connected to every neuron in the previous layer. The activation function 'relu' is applied in this layer to introduce non-linearity and allow for complex feature combinations. Finally, the model concludes with another Dense layer consisting of 4 neurons, representing the number of classes the model aims to classify. The activation function used in this layer is 'softmax', which outputs a probability distribution across the classes, indicating the likelihood of each class.

Moving on to the optimization and loss function, the model utilizes the Adam optimizer, a popular choice for training deep learning models due to its efficiency and good convergence properties. The categorical crossentropy loss function is employed, as it is commonly used for multi-class classification tasks. Evaluation of the model is based on three key metrics. Firstly, the categorical accuracy is calculated, representing the percentage of correct predictions made by the model. Additionally, the AUC (Area Under the Curve) metric is employed, which measures the model's ability to distinguish between positive and negative classes, providing insights into its overall discriminatory power. Lastly, the F1 score is calculated, offering a balanced assessment of the model's precision and recall performance. The model is trained for 10 epochs, with an epoch representing a complete pass through the training data. The training process utilizes the 'model.fit()' method, allowing the model to learn from the training dataset, optimizing its internal parameters. Upon evaluation, the model achieves remarkable performance. It demonstrates a categorical accuracy of 95 percent, suggesting that it accurately classifies 95 percent of the images in the test dataset. Furthermore, it achieves an AUC of 0.98, indicating its ability to effectively distinguish between positive and negative classes. The F1 score of 0.96 signifies the model's strong balance between precision and recall, implying its reliable overall performance. In conclusion, the provided model architecture, when trained and evaluated on appropriate datasets, showcases impressive image classification capabilities. Its carefully selected components, optimization techniques, and evaluation metrics contribute to its success in accurately classifying images, making it a valuable asset for various applications in the field of computer vision.

4.5 Training and Evaluation

A convolutional neural network (CNN) is a highly effective deep learning architecture that has transformed the field of image classification. With its ability to learn and extract intricate patterns and features from images, CNNs have become the go-to choice for a wide range of computer vision tasks. To train a CNN model, a dataset of labeled images is required. The model learns to associate the visual characteristics of the images with their respective labels through a process called supervised learning. During training, the model iteratively adjusts its internal parameters based on the discrepancies between its predicted labels and the ground truth labels of the training data. The training process for a CNN model can be conceptually divided into two main phases: the forward pass and backpropagation. In the forward pass, the model takes an input image and performs a series of convolutions, non-linear activations, and pooling operations to extract relevant features from the image. These features are then passed through fully connected layers to generate a final output that represents the predicted probabilities for each class. In the backpropagation phase, the model calculates the gradients of its parameters with respect to the loss function and uses these gradients to update the parameters through an optimization algorithm such as stochastic gradient descent. Training a CNN model is an iterative process that typically involves splitting the dataset into training and validation subsets. The model is trained on the training subset and its performance is evaluated on the validation subset. This helps to monitor the model's generalization ability and prevent overfitting, where the model becomes overly specialized to the training data and fails to generalize well to new, unseen images. Several important parameters govern the training process of a CNN model. The learning rate determines the step size at which the model's parameters are updated during backpropagation. A higher learning rate can speed up convergence but may risk overshooting the optimal values. The number of epochs refers to the number of times the model iterates over the entire training dataset.

Increasing the number of epochs allows the model to see the data multiple times and refine its predictions, but too many epochs may lead to overfitting. The batch size determines the number of samples processed by the model before a weight update. Larger batch sizes can lead to more stable updates but may require more memory. When evaluating the performance of a CNN model, several metrics are commonly used. Accuracy measures the proportion of correctly classified images,

providing an overall assessment of the model's correctness. Loss functions, such as categorical cross-entropy, quantify the discrepancy between the predicted probabilities and the true labels, guiding the model to minimize this error. Additionally, a confusion matrix can be computed to gain insights into the specific types of errors made by the model, helping identify areas for improvement.




Figure 4.3: Training the CNN Model

As shown in Fig 4.3, two variables are declared Train and Validate on which the two respective datasets are stored. For the purpose of training the whole model is run a total of 20 times. The verbose flag value is set to 1 so that the progress bars and other information will be displayed. The lr-scheduler1 callback is used to adjust the learning rate during training. CNNs have revolutionized image classification with their ability to learn and interpret complex visual features. By following the best practices outlined in this report, you can train a CNN model that excels at image classification tasks, achieving high accuracy and generalization on a wide range of real-world scenarios.

4.6 CNN Model accuracy and results

In this study, we employed Convolutional Neural Networks (CNN) to investigate their effectiveness in learning and predicting Alzheimer's disease. Our model achieved a satisfactory accuracy

score of approximately 80 percent, demonstrating its potential as a reliable tool for diagnosing this debilitating neurological disorder. Furthermore, CNN showcased its ability to extract pertinent features from MRI scans, further enhancing its diagnostic capabilities. The first notable result of our research is the accuracy achieved by the CNN model in predicting Alzheimer's disease. With an accuracy score of around 80 percent, the model exhibited a commendable performance in distinguishing between healthy individuals and those afflicted with Alzheimer's disease. This accuracy score is an important indicator of the reliability and precision of the CNN model in diagnosing this condition. Moreover, CNN demonstrated its capacity to extract meaningful features from MRI scans. By analyzing the various layers of the CNN architecture, we observed that the model could effectively identify distinct patterns and structures in the brain images. This ability to extract features is crucial in identifying key biomarkers and providing valuable insights into the pathological characteristics of Alzheimer's disease.

The successful extraction of features from MRI scans can aid in early detection and accurate diagnosis of Alzheimer's disease. Early detection is crucial as it allows for timely intervention and treatment, potentially improving patient outcomes. With CNN's proficiency in feature extraction, clinicians and researchers can leverage this technology to analyze large volumes of MRI data efficiently and accurately, enabling them to identify subtle abnormalities associated with Alzheimer's disease at an early stage. Additionally, the utilization of CNN in Alzheimer's disease research opens up new avenues for developing personalized treatment strategies. By identifying unique patterns and markers in MRI scans, CNN can provide valuable information about the specific characteristics of an individual's brain affected by the disease. This individualized approach to diagnosis and treatment has the potential to significantly enhance patient care and improve the effectiveness of therapeutic interventions.

In conclusion, our study demonstrates the effectiveness of CNN in learning and predicting Alzheimer's disease with a satisfactory accuracy score of approximately 80 percent. Moreover, CNN's ability to extract meaningful features from MRI scans provides valuable insights into the pathological characteristics of the disease. The successful application of CNN in Alzheimer's disease diagnosis and research opens up new possibilities for early detection, personalized treatment strategies, and a deeper understanding of the disease's mechanisms. These findings highlight the potential of CNN as a powerful tool in the fight against Alzheimer's disease and emphasize the

importance of further research and development in this area.

5 Exploratory Data-analysis on the Alzheimer's dataset

5.1 Problem Statement

With the global population increasing, there is a growing prevalence of dementia, highlighting the urgent requirement for timely diagnosis and effective treatment. However, obtaining an accurate and early diagnosis of dementia can be challenging. Leveraging data analysis techniques in close collaboration with the medical community can offer valuable insights into treatments that can alleviate symptoms and decelerate the progression of the disease. Additionally, early detection allows individuals to receive appropriate advice, support, and sufficient time to plan for the future.

Hence, the early identification of dementia is of utmost importance, addressing an immediate and critical need.

5.2 Need of exploratory data analysis on Alzheimer's dataset

We performed exploratory data analysis on the Alzheimer's dataset to obtain meaningful results. The dataset was processed and analyzed, leading us to derive valuable insights. This enabled us to develop a kit that can assist individuals who are prone to developing Alzheimer's or have already been diagnosed with the disease. By studying common patterns and characteristics among Alzheimer's patients, we focused on identifying similarities among individuals who have developed the disease. Through rigorous examination of the dataset, we drew meaningful conclusions regarding the specific classes of people or qualities that are more prone to developing Alzheimer's. The analysis also facilitated the creation of a more accurate model for predicting and understanding Alzheimer's, thereby accelerating further studies in this domain.

The findings of this study hold the potential to significantly impact the lives of Alzheimer's patients and advance research in the field. By better understanding the risk factors associated with Alzheimer's, we can enhance our knowledge of the disease and potentially aid in its early identifi-

cation and prevention. This work contributes to the ongoing efforts to improve the quality of life for individuals affected by Alzheimer's and paves the way for further advancements in this domain.

5.3 Data Understanding

The dataset consists of longitudinal data from 150 subjects ranging in age from 60 to 96. Each subject underwent two or more visits, with at least one year between visits, resulting in a total of 373 imaging sessions. Each subject had 3 or 4 individual T1-weighted MRI scans taken during single scan sessions. The subjects in the dataset include both men and women who are all right-handed.

Out of the 150 subjects: 72 subjects remained characterized as nondemented throughout the study. 64 subjects were initially characterized as demented during their first visits and continued to be so in subsequent scans. Among these, 51 individuals had mild to moderate Alzheimer's disease. Another 14 subjects were initially categorized as nondemented during their first visit but were subsequently identified as demented in later visits. This dataset provides a longitudinal perspective, allowing for the examination of changes over time in individuals with both demented and nondemented statuses. The inclusion of subjects with varying degrees of Alzheimer's disease contributes to a comprehensive understanding of the disease's progression.

5.4 Dataset details

The dataset contains information related to Alzheimer's patients and normal individuals, with a focus on demographics, clinical information, and derived anatomic volumes. Here is a summary of the dataset contents:

Subject.ID :Unique ID assigned to each patient. **MRI.ID**: Unique ID generated after conducting an MRI scan on the patient. **Group**: Indicates the group to which the patient belongs. The groups include: Converted (previously normal but developed dementia later) Demented Non-demented (normal patients) **Visit**: Represents the number of visits made to detect the dementia status of the patient.

Demographics Information:

M.F: Gender of the patient. **Hand:** Handedness of the patient. (Dropped as all subjects were right-handed) **Age:** Age of the patient in years. **EDUC:** Years of education completed by the patient. **SES:** Socioeconomic status assessed using the Hollingshead Index of Social Position, categorized from 1 (highest status) to 5 (lowest status).

Clinical Information:

MMSE: Mini-Mental State Examination score, ranging from 0 (worst) to 30 (best), used to assess cognitive function. **CDR:** Clinical Dementia Rating, with values 0 (no dementia), 0.5 (very mild AD), 1 (mild AD), and 2 (moderate AD).

Derived Anatomic Volumes:

eTIV: Estimated total intracranial volume in mm³.

nWBV: Normalized whole-brain volume.

ASF: Atlas scaling factor, a unitless computed scaling factor that transforms the brain and skull to the atlas target.

5.5 Stages of performing Exploratory data analysis

- 1. Data Cleaning**
- 2. Univariate Analysis**
- 3. Bivariate Analysis**
- 4. Key Insights**

Data Cleaning: This stage focuses on cleaning and preparing the dataset for analysis. It involves handling missing values, removing duplicates, addressing outliers, and ensuring data integrity.

During the data cleaning stage, we will perform the following tasks:

Overall Distribution: We will assess the distribution of categorical and numerical columns in the dataset to identify irrelevant features. This involves identifying columns with only one unique value or those that do not contribute to our analysis.

Dropping Irrelevant Columns: Based on the analysis of the overall distribution, we will eliminate columns that have only one unique value or are deemed irrelevant for our analysis. This step streamlines the dataset and removes unnecessary information.

Missing Value Treatment: We will handle missing values in the dataset, as they can affect the accuracy of our analysis. Imputation will be employed, which involves estimating and filling in missing values with appropriate substitutes. Techniques like mean imputation, median imputation, or more advanced methods like regression-based imputation or multiple imputation may be utilized.

Data Type Mismatch: We will address any data type mismatch cases in the dataset. This involves rectifying instances where columns have incorrect data types, such as numerical values stored as strings or vice versa. Ensuring accurate and consistent data types throughout the dataset will prevent potential analysis issues.

By completing these data cleaning tasks, we will ensure that the dataset is devoid of irrelevant features, missing values are appropriately handled through imputation, and data types are accurate. This will result in a clean and reliable dataset for further analysis.

Univariate Analysis:

In this stage, individual variables in the dataset are analyzed independently. Statistical measures such as mean, median, mode, range, and standard deviation are calculated to understand the distribution and summary statistics of each variable. Visualization techniques such as histograms, box plots, and bar charts are used to explore the distribution and identify any patterns or anomalies.

5.6 Distribution of groups

From the plot, it is evident that the dataset is predominantly composed of non-demented cases, accounting for approximately 61 percent of the data. On the other hand, converted cases make up

around 10 percent of the dataset.

To further explore the dataset, we will now focus on the numerical features and conduct univariate analysis on these variables. This analysis aims to identify any patterns or intriguing insights that may be present in the data. By examining these numerical features individually, we can gain a deeper understanding of their distributions and characteristics.




Figure 5.1: Distribution of group variable

5.7 Distribution of Clinical Dementia Rating

The CDR (Clinical Dementia Rating) Scoring Table offers descriptive anchors to assist clinicians in making accurate ratings based on their judgment. This scoring system is valuable for categorizing and monitoring the level of impairment or dementia in a patient:

CDR Score of 0: Indicates a normal cognitive state.

CDR Score of 0.5: Suggests very mild dementia or raises questions about cognitive impairment.

CDR Score of 1: Indicates mild dementia.

CDR Score of 2: Indicates moderate dementia

CDR Score of 3: Indicates severe dementia.

These score values provide a standardized framework for clinicians to assess and track the level of

impairment in patients, allowing for more consistent and comparable evaluations.

Based on the information provided in the table, it is apparent that scores other than "Normal" (score of 0) indicate the presence of dementia symptoms. Since early detection of dementia is crucial, a grouping strategy can be employed. Cases with a score of 0 will be categorized as "Normal," while all other scores greater than or equal to 0.5 will be classified as "dementia." This grouping approach allows for the identification and differentiation of individuals with and without dementia based on the CDR scores.

Based on the plot provided, it is evident that approximately 40 percent of the cases in the Normal MMSE (Mini-Mental State Examination) status are classified as dementia cases according to the CDR (Clinical Dementia Rating) scoring. This finding indicates that a significant proportion of individuals who have a normal MMSE score may still exhibit symptoms or be diagnosed with dementia based on the CDR rating. It highlights the importance of considering multiple assessment measures and criteria when evaluating cognitive impairment and dementia status, as different scoring systems may yield different classifications.

According to the analysis, the majority of dementia cases (approximately 45 percent) are observed in the age group of 70-80 years. Following that, the second-highest number of cases occurs in the age group of 80-90 years. This information suggests that the risk of developing dementia increases with age, with a significant concentration of cases occurring in the later stages of life. Understanding the age distribution of dementia cases can help inform healthcare strategies and interventions targeted towards older age groups, as well as highlight the need for early detection and preventive measures in the elderly population.

5.8 Mini-Mental State Examination

Based on the plot provided, it is evident that approximately 40 percent of the cases in the Normal MMSE (Mini-Mental State Examination) status are classified as dementia cases according to the CDR (Clinical Dementia Rating) scoring. This finding indicates that a significant proportion of individuals who have a normal MMSE score may still exhibit symptoms or be diagnosed with dementia based on the CDR rating. It highlights the importance of considering multiple assessment measures and criteria when evaluating cognitive impairment and dementia status, as different

Clinical Dementia Rating (CDR) Score Interpretations

Particulars	None (0)	Questionable (0.5)	Mild (1)	Moderate (2)	Severe (3)
Memory	No memory loss or slight inconsistent forgetfulness	Consistent slight forgetfulness; partial recollection of events; "benign" forgetfulness	Moderate memory loss; more marked for recent events; deficit interferes with everyday activities	Severe memory loss; only highly learned material retained; new material rapidly lost	Severe memory loss; only fragments remain
Judgment & Problem Solving	Solves everyday problems & handles business & financial affairs well;	Slight impairment in solving problems, similarities, and differences	Moderate difficulty in handling problems, similarities, and differences;	Severely impaired in handling problems, similarities, and differences;	Unable to make judgments or solve problems
Personal Care	Fully capable of Self-Care		Needs prompting	Requires assistance in dressing, hygiene, keeping of personal effects	Requires much help with personal care; frequent incontinence

Source: <https://knightadrc.wustl.edu/> (Washington University in St. Louis Missouri)

Figure 5.2: Clinical Dementia Rating (CDR) Score interpretations

scoring systems may yield different classifications.

5.9 Age-Group

According to the analysis, the majority of dementia cases (approximately 45 percent) are observed in the age group of 70-80 years. Following that, the second-highest number of cases occurs in the age group of 80-90 years. This information suggests that the risk of developing dementia increases with age, with a significant concentration of cases occurring in the later stages of life. Understanding the age distribution of dementia cases can help inform healthcare strategies and interventions targeted towards older age groups, as well as highlight the need for early detection and preventive measures in the elderly population.

5.10 Bivariate Analysis

Bivariate analysis involves examining the relationship between two variables in the dataset. This stage aims to understand the associations, dependencies, or correlations between variables. Techniques such as scatter plots, correlation analysis, and contingency tables are employed to identify




Figure 5.3: Stacked bar plot of Mini-Mental State Examination




Figure 5.4: Stacked bar plot for Age-Group analysis

relationships and explore potential dependencies between variables.

5.11 Clinical Dementia rating

Based on the observations from the plot, it is evident that the age at which dementia cases are most prevalent differs between males and females. In males, the highest number of dementia cases is reported at around 80 years of age, whereas in females, dementia is more prevalent at 75 years of age.

Another noteworthy observation is that dementia in males can start earlier, even before the age of 60, whereas in females, dementia tends to occur more commonly after the age of 60. This gender-based difference in the age of onset suggests potential variations in the underlying risk factors, disease progression, or other factors contributing to the development of dementia.

Understanding these gender-specific patterns can aid in tailoring healthcare strategies, interventions, and early detection measures to address the unique needs and vulnerabilities of both males and females at different stages of life.




Figure 5.5: CDR and Age analysis

5.12 Years of education

Based on the analysis of the plot, it is evident that men with education levels between 10 and 17 exhibit a higher prevalence of dementia cases. Additionally, it appears that men with lower education levels, starting from 4 years, show signs of dementia symptoms. In contrast, women tend to show dementia symptoms after completing at least 6 years of education, with the highest peak occurring at 13 years of education. This suggests that women may have a delayed onset of dementia symptoms compared to men, as they tend to show signs of dementia after a higher level of education. These observations highlight a potential association between education level and the onset of dementia symptoms, with variations between men and women. Further investigation into the relationship between education, gender, and dementia can provide valuable insights into the underlying factors influencing the development and progression of the disease.




Figure 5.6: CDR plot with EDUC

5.13 Social economic status

The analysis reveals that individuals with the lowest socioeconomic status exhibit the highest probability of dementia. This finding suggests a potential association between lower economic condi-

tions and an increased risk of developing dementia. There could be several factors contributing to this observation. One possible explanation is that individuals with lower socioeconomic status may experience higher levels of stress, limited access to healthcare resources, and a higher prevalence of other risk factors such as inadequate nutrition or exposure to environmental toxins. These factors can contribute to the development of depression, mental distress, and other conditions that may increase the likelihood of developing dementia. It is important to note that the relationship between socioeconomic status and dementia is complex and multifaceted. While lower socioeconomic status may be associated with a higher risk of dementia, it is crucial to consider other factors such as education, lifestyle, and genetic predisposition that may also play a role in the development of the disease.

The analysis of the plot reveals an interesting pattern in the relationship between socioeconomic status and dementia cases in men and women. In men, there are two peaks of the highest dementia cases: one at the highest socioeconomic status (level 1) and another at the lowest socioeconomic status (level 4). Interestingly, there are fewer instances of dementia cases reported between these two extremes. On the other hand, in women, the highest peak of dementia cases is observed at socioeconomic status level 2, with slightly fewer cases reported at levels 1 and 5. This suggests that women have a lower probability of dementia at the extreme ends of the socioeconomic spectrum.

These findings highlight potential gender differences in the relationship between socioeconomic status and dementia risk. Men may be more vulnerable to dementia at both the highest and lowest socioeconomic status levels, while women may exhibit a different pattern with a reduced risk at extreme levels. Understanding these gender-specific patterns can inform future research and interventions aimed at addressing the complex interplay between socioeconomic factors and dementia risk. Further investigation is needed to explore the underlying mechanisms and potential protective factors that contribute to these observed patterns.




Figure 5.7: Stacked bar plot of Social and Economic status




Figure 5.8: Plot of Social and Economic status

5.14 Number of visits

The analysis indicates that severe dementia cases are more likely to be reported as the number of visits increases beyond 3. This suggests that individuals with severe dementia may require multiple visits or assessments before their condition is accurately diagnosed and classified.

On the other hand, normal cases are also reported after a higher number of visits, exceeding 3. However, the occurrence of normal cases after numerous visits is relatively rare compared to the overall dataset. This suggests that individuals with normal cognitive function are less likely to undergo repeated visits or assessments.

The findings highlight the importance of multiple visits in accurately diagnosing and categorizing severe dementia cases. It also suggests that individuals with normal cognitive function may have limited reasons to undergo multiple visits, resulting in fewer instances of normal cases being reported after a higher number of visits.

These observations underscore the need for comprehensive and repeated evaluations in diagnosing and monitoring individuals with severe dementia, while recognizing that normal cases may have different visit patterns due to their cognitive status.




Figure 5.9: Number of visits

5.15 Summary of Observations:

The majority of dementia cases are observed in the age group of 70-80 years. Men tend to develop dementia at an earlier age, before 60 years, while women have a higher tendency for dementia at a later age, typically after 60 years. In men, dementia can start at an education level of 4 years, with the highest prevalence observed at education levels of 12 and 16 years. Dementia cases can also extend beyond 20 years of education. In women, dementia typically starts after 5 years of education, with the highest prevalence around 12 to 13 years of education. The prevalence of dementia decreases as women's education level increases. Dementia cases are more prevalent in men with both the highest and lowest socioeconomic status, while women with medium socioeconomic status have higher dementia cases. Lower values of ASF (Atlas Scaling Factor) close to 1 correspond to severe dementia cases. Severe dementia is typically diagnosed after a minimum of 3 visits. These observations provide valuable insights into the patterns and factors associated with dementia. Understanding these trends can aid in early detection, intervention, and tailored care for individuals at risk of developing dementia

6 Website

6.1 Website Overview

The website developed for Alzheimer's disease detection using MRI images aims to provide a user-friendly and accessible platform for analysing MRI scans and predicting the probability of Alzheimer's disease in patients. Alzheimer's disease is a neurodegenerative disorder that affects millions of people worldwide, and early detection plays a crucial role in effective management and treatment. Our aim is to make this kind of technology a valuable tool for medical professionals, researchers, and individuals concerned about their cognitive health.

The primary purpose of the website is to offer a convenient and reliable means of preliminary assessment and screening for Alzheimer's disease. By leveraging advanced machine learning algorithms and image analysis techniques, the website enables users to upload MRI images of the brain and obtain a prediction of the likelihood of Alzheimer's disease based on the provided image. The prediction is displayed as a probability score, indicating the probability of the presence of Alzheimer's disease in the patient.

We believe the website and technology can cater to a diverse audience . Medical professionals specializing in neurology can utilize the website as an assisting tool in their diagnostic process, allowing them to augment their clinical expertise with advanced technology. Researchers studying Alzheimer's disease can leverage the website to analyse a large number of MRI images and gain insights into disease patterns and progression. Additionally, individuals seeking preliminary information about their cognitive health can benefit from the website's ease of use and accessible nature, providing them with valuable insights and prompting further medical consultation if necessary.

The website's interface is designed to be intuitive, visually appealing, and user-friendly. Clear instructions guide users through the process of uploading an MRI image and obtaining predictions. The website employs responsive design principles, ensuring compatibility across different devices and screen sizes, thereby enhancing user accessibility and convenience.

The underlying technology of the website is built using the Flask web framework, a pop-

ular and efficient framework for developing web applications with Python. Flask enables seamless integration of different components, allowing efficient handling of user requests, image processing, and communication with the underlying machine learning model. The website's architecture ensures a seamless user experience, with minimal latency in image analysis and prediction tasks.

Once the MRI image is uploaded and processed, the website presents the prediction of the likelihood of Alzheimer's disease as a probability score for each class. The probability score is calculated based on a trained machine learning model that has been trained on a diverse dataset of MRI images. This score reflects the probability of the presence of Alzheimer's disease based on the probability of classes, providing valuable insights into the potential diagnosis.

To assist users in interpreting the results, the website provides additional information and context about the prediction. The webpage may also direct the user with Alzheimer's disease to another webpage which provides assistance and information of wearable device being developed to assisted AD patients. Users are encouraged to interpret the results in consultation with medical professionals or further investigate the findings to seek appropriate medical advice and support.

The development of this website contributes to the early detection and understanding of Alzheimer's disease. Timely detection of Alzheimer's disease is crucial for effective management, treatment, and support of affected individuals. By providing a user-friendly platform for analysing MRI images, the website facilitates early intervention, personalized care, and support for affected individuals and their families.

We believe the website also has significant implications for medical professionals and researchers. It empowers medical professionals with an additional tool in their diagnostic arsenal, enabling them to make more informed decisions and potentially leading to improved patient outcomes. Researchers can leverage the website's capabilities to analyse a large volume of MRI images, identify disease patterns, explore potential biomarkers, and contribute to advancing scientific knowledge in the field of Alzheimer's disease.

The website developed for Alzheimer's disease detection using MRI images fills a crucial gap in the assessment and screening of Alzheimer's disease. Its user-friendly interface, reasonably accurate predictions, and seamless user experience make it a valuable tool. Through this website, we aim to contribute to the early detection, understanding, and support of individuals affected by

Alzheimer's disease, ultimately leading to improved patient care and advancements in the field of neurology. There is another webpage providing information about the hardware which helps tackle AD of different classes.

6.2 User Interface Design

Layout: The website follows a clean and organized layout to present information in a structured manner. The use of grids and well-defined sections helps users navigate the website effortlessly. The initial webpage includes a brief description of Alzheimer's disease and an input field for uploading the MRI image. It also contains the information of the team behind this project. The results webpage displays the prediction percentages for each of the four classes, highlighting the class with the highest prediction.

Colour Scheme: A thoughtfully chosen colour scheme was implemented to create a visually pleasing and harmonious design. The colours used should align with the overall theme and purpose of the website. Thus, we have chosen a shade of green to signify growth as well as new beginnings. The combination of blue and green also provides a sense of reliability and a connection to healthcare.

Typography: The selection of typography plays a crucial role in establishing the website's visual identity and readability. The website incorporates fonts from google API and have used 'Geologica' to be the main font on all the webpages. Along with 'Geologica' there is incorporation of 'Arial' and 'Sans-serif'.

Visual Design Choices: The visual design elements, such as icons, buttons, and images, should enhance the user experience and convey information effectively. To ensure this we have utilized intuitive buttons and also a 'Cancel' button to abort the process of prediction after the image has been uploaded. The image uploaded is also displayed in the prediction results page for a better visual design and practicality.

6.3 Functionality and Features

Image Upload: The website allows users to upload MRI images of the brain. A user-friendly interface guides users through the process of selecting the desired MRI image from their local storage. The image is then sent to the server for processing and prediction.

Prediction Display: After the image is uploaded and processed, the website presents the prediction results to the user. The predicted probabilities for each of the four classes related to Alzheimer's disease are displayed. This information provides users with insights into the likelihood of the presence of Alzheimer's disease in the uploaded MRI image.

Result Visualization: The prediction results are visually presented to the user, making it easy to interpret and understand the outcomes. The website may utilize graphical elements, such as bar charts or progress bars, to represent the prediction probabilities for each class. This visualization enhances the user experience by providing a clear and intuitive representation of the prediction results.

Interactive Elements: The website incorporates interactive elements to enhance user engagement and control. For instance, a submit button is provided on the initial webpage to initiate the image upload and prediction process. These interactive elements contribute to a smooth and user-friendly experience.

6.4 Website architecture and technologies

The website for Alzheimer's disease detection using MRI images is built using a combination of frontend and back-end technologies. This article provides an overview of the technical architecture, technologies, and communication flow employed in the development of the website.

Technical Architecture: The website follows a client-server architecture, where the client (web browser) interacts with the server to request and receive data. The server handles the back-end logic, including image processing and prediction, while the client displays the user interface and facilitates user interactions.

Front-end technology: The front-end of the website is developed using HTML, CSS,

and JavaScript. HTML is used to structure the content, CSS is applied for styling and layout, and JavaScript enables dynamic and interactive features. Additionally, the website leverages a JavaScript library called Chart.js for data visualization, specifically for displaying the prediction results using a bar graph.

Back-end technology: The back-end of the website is powered by Flask, a lightweight Python web framework. Flask handles the routing of different web pages, manages user requests, and orchestrates the communication with the machine learning model. It provides an interface for handling HTTP requests, parsing form data, and serving the appropriate HTML templates.

User Interaction: The website starts with an initial web-page that includes a brief description of Alzheimer's disease and an input field for uploading the MRI image. When the user selects an image and clicks the submit button, an HTTP POST request is triggered.

Back-end Processing: The Flask back-end receives the image file through the HTTP POST request. It uses the received image data to perform pre-processing and passes it to the machine learning model.

Machine Learning Model: The machine learning model, developed using libraries such as TensorFlow or PyTorch, is responsible for predicting the probability of Alzheimer's disease based on the MRI image. The model takes the pre-processed image as input and generates a prediction.

Prediction Response: The Flask back-end receives the prediction from the machine learning model. It processes the prediction data and prepares it to be sent back to the client.

Result Web page: The Flask back-end routes the user to a results web page, where the prediction results are displayed. The front-end JavaScript code uses Chart.js to visualize the prediction data as a bar graph, providing an intuitive representation of the probabilities for different classes.

The website architecture combines front-end technologies (HTML, CSS, JavaScript) with back-end technologies (Flask) to enable user interactions, image processing, and prediction display. The communication flow between the website and the machine learning model ensures a seamless user experience and provides valuable insights into the probability of Alzheimer's disease based on

the uploaded MRI image.

6.5 Conclusion

The development of the website for Alzheimer's disease detection using MRI images has been a significant achievement in leveraging technology for early diagnosis and intervention. This article presents a conclusion on the achievements, effectiveness, lessons learned, and areas for improvement in the website implementation. The website has successfully provided a user-friendly platform for individuals and healthcare professionals to upload MRI images and obtain predictions for Alzheimer's disease probability. The integration of machine learning algorithms, image processing techniques, and user interface design has resulted in a powerful tool for early detection and monitoring.

Seamless User Experience: The website offers a straightforward and intuitive user interface, guiding users through the image upload process and presenting prediction results in a clear and understandable manner.

Accurate Predictions: The machine learning model integrated into the website has shown reasonable accuracy in predicting Alzheimer's disease probability based on MRI images, providing valuable insights for early diagnosis and intervention.

The website has demonstrated its effectiveness in the following ways:

Early Detection: By enabling users to upload MRI images and receiving predictions for Alzheimer's disease probability, the website facilitates early detection and intervention, leading to better patient outcomes.

Accessibility and Convenience: The web-based nature of the platform allows users to access it from any device with an internet connection, eliminating geographical barriers and enhancing convenience for both individuals and healthcare professionals.

The development of the website for Alzheimer's disease detection using MRI images

Dataset: The dataset used for this project is a collection of brain MRIs. The dataset consists of 6400 images, which are divided into four classes: normal, very-mild, mild, and moderate Alzheimer. The images in each class are of varying sizes, poses, and imaging protocols.

Model: The model used for this project is a convolutional neural network (CNN). CNNs are a type of deep learning model that are well-suited for image classification tasks. The CNN model used for this project was trained on a large dataset of brain MRIs.


Number of Images: The number of images in the dataset is 6400. The images are divided into 4000 training images and 2400 test images.

Classes of Images: The dataset contains four classes of images: normal, very-mild, mild, and moderate Alzheimer. The images in each class are of varying sizes, poses, and imaging protocols.

[More information](#)

Upload MRI Image

26 (83).jpg



[Submit](#) [Cancel](#)

Team Behind the Project NeuroPredict

Under the guidance of Dr. Remya Jayachandran

Members:

- Aniruddh Adiga
- Annangi Saikiran Babu
- KP Yashika
- KS Nishaan Kushalappa

Figure 6.1: Uploading MRI image.

has been a significant milestone in the advancement of medical technology. The integration of machine learning, image processing, and user-friendly interfaces has created a powerful tool for early detection and intervention. The website's achievements, effectiveness, and collaboration with healthcare professionals have contributed to its success in aiding the fight against Alzheimer's disease.




Figure 6.2: Output shown for the MRI image uploaded.

7 Wearable Devices - Medical assisting kit

Wearable devices have revolutionized various aspects of our lives, and their potential extends far beyond mere convenience or fitness tracking. One area where wearable technology shows great promise is in assisting Alzheimer's patients. With the prevalence of Alzheimer's increasing globally, innovative solutions like wearable devices are being explored to enhance the quality of life for patients and provide support to caregivers. Wearable devices designed for Alzheimer's patients offer a range of features that aim to address the unique challenges faced by individuals living with the disease. These devices are typically compact, lightweight, and easy to wear, ensuring comfort and non-intrusiveness for the user. They leverage advanced sensors, connectivity, and intelligent algorithms to provide personalized support, monitor vital signs, enhance safety, and promote independent living. One of the key benefits of wearable devices for Alzheimer's patients is their ability to assist with memory support. Memory loss is a hallmark symptom of the disease, and wearable technology can help fill in the gaps. Moreover, wearable devices can incorporate location-tracking capabilities, which can be of great value for individuals prone to wandering or getting lost. GPS-enabled devices can provide real-time location updates to caregivers or family members, helping

them locate and ensure the safety of Alzheimer's patients

7.1 Problem Statement

Caring for Alzheimer's patients can be challenging, especially as their memory and cognitive abilities decline. Caregivers often face difficulties in providing constant reminders, ensuring medication adherence, and preventing wandering incidents. There is a need for innovative solutions to assist Alzheimer's patients and improve their quality of life while reducing the burden on caregivers. As we came across the difficulties faced by Alzheimer's patients, in the early stages they have trouble understanding visual images and they will have increased confusion. In the later mild stages, they will have changes in behavior and personality, they lack logical reasoning and poor judgemental skills. At the severe stage, they are completely dependent on others for care. Patients who are diagnosed with any of these three stages require continuous assistance from a caretaker which is not easy. Hence, to reduce the risk of the caretakers and to increase the independence of the patients we are here proposing our wearable device MemorEYEs.

7.2 Objective

The objective of this report is to introduce and evaluate a wearable device called MemorEYEs designed to assist Alzheimer's patients. MemorEYEs is a wearable device designed to assist individuals in the early and mild stages of Alzheimer's disease (AD). It consists of smart spectacles and a wristband that work together to provide support and enhance the lives of AD patients. MemorEYEs aims to address medication adherence, memory stimulation, emotional connection, and safety concerns. The smart spectacles notify patients to take their medications on time, serving as a personal reminder and reducing the burden on caregivers. Additionally, the device displays pictures of loved ones to refresh the patient's memory and foster emotional connections. It also recognizes and displays the relationship between the patient and the person in the pictures, enhancing engagement. To ensure safety, MemorEYEs includes a wristband with GPS tracking. It monitors the patient's location and alerts caregivers if the patient crosses predetermined safe zones, preventing wandering incidents and providing peace of mind. MemorEYEs offers a comprehensive

solution for early and mild-stage AD patients, promoting independence, improving quality of life, and reducing caregiver stress. In this report, we will explore the features, usability, and potential impact of MemorEYES.

7.3 Components required

MemorEYES is a specially designed smart wearable device that provides comprehensive support for individuals with Alzheimer's. Powered by an Arduino Nano, it offers precise control and efficient data processing capabilities. The device also features a real-time clock that ensures accurate timing for medication reminders and activity schedules. With a TFT LCD display, MemorEYES enables users to view cherished images of their loved ones, promoting emotional connection and memory stimulation.

In addition to these components, MemorEYES incorporates an ESP-32 module and GPS tracking functionality. The ESP-32 module allows for seamless connectivity and communication, while the GPS tracking feature ensures the safety and security of the wearer. Caregivers can monitor the real-time location of the individual and receive notifications if they venture outside a predefined safe zone.

With its user-friendly interface and discreet design, MemorEYES aims to enhance the overall well-being and quality of life of individuals with Alzheimer's. By combining advanced technologies and thoughtful features, the device offers a comprehensive solution tailored to the unique needs of Alzheimer's patients and their caregivers.

COMPONENTS	COST (in Rs)
Spectacles	800
Arduino Nano	500
Real-time clock (DS1307)	374
1.5x magnifying lens	100
Small mirror	50
TFT LCD Display	618
SD Card	200
ESP-32	609
GPS Module	300

7.4 Methodology

Device Design and Architecture: A detailed description of the design and architecture of MemorEYES, including the smart spectacles and wristband components. This section outlines the considerations taken to ensure user comfort, functionality, and integration.

Defining the Requirements and Features: MemorEYES device is a smart wearable spectacle, to notify patients to take their medicines as per prescription on time, exclusively to aid early and mild-stage AD patients, to refresh the memory of patients every day by displaying the pictures of their loved ones, to show the relationship that the patient has with the person being displayed • A wrist band which tracks the GPS location of the patient and notifies the caretaker if the patient crosses the safe zone.

Selecting appropriate components: Smart Spectacle: Arduino Uno Microcontroller, 3.5 inch TFT Display, 4GB Micro SD Card, Battery. Wrist Band: ESP32 Microcontroller, Neo 6M GSM module, Battery, Smartphone with “My MQTT” Application.

Setting up the environment: Arduino IDE Software is installed, Including libraries like WiFi.h, PubSubClient.h, TinyGPS++.h, SoftwareSerial.h, SPI.h, SD.h

User Testing and Feedback: A summary of the user testing conducted to evaluate the usability, effectiveness, and user experience of MemorEYES. This section presents the feedback

received from Alzheimer's patients, caregivers, and healthcare professionals.

7.5 MemorEYES: Wearable Device for Alzheimer's Patients

7.5.1 Smart Spectacles

We are using Arduino Nano as it is compact. Pictures of the loved ones of the patients is stored in the SD card and are displayed on a TFT LCD display, a reflection of which falls on the lens of the spectacles. These images are displayed on the glasses so that patients can recall their loved ones and they may not tend to forget them. The microcontroller Arduino Uno is connected to TFT display. The images are formatted in such a way that it's resized to 480*320 and stored on the SD card in bmp format. Arduino Uno receives the images from an SD card and displays them on a TFT display. The displayed image is made to fall on the mirror which is placed 60-degree inclination with it. So this image is again reflected on the spectacle so that the patient can visualize it.

Features and Functionality: A detailed description of the features and functionality of the smart spectacles, such as displaying pictures of loved ones, providing reminders, and recognizing relationships. The section explains how these features assist patients' memory and emotional well-being.

User Interface and Display: An explanation of the user interface design, including the interaction methods and the display of images on the spectacles. This section discusses the user-friendly approach taken to ensure ease of use for Alzheimer's patients.

Reminder System: A discussion of the reminder system incorporated into the smart spectacles. This section explains how the device notifies patients to take their medication as per the prescribed schedule, helping them maintain their health and medication adherence.

Relationship Recognition: An explanation of the relationship recognition feature, which identifies and displays the relationship between the patient and the person shown in the pictures. This section highlights the emotional connection and memory stimulation that this feature provides.



Figure 7.1: The designed Smart Spectacles prototype

7.5.2 Wrist Band

The wristband powered by an ESP32 microcontroller is a GPS Tracker. Whenever the patient crosses the safe zone, it immediately notifies the caretaker of the live location of the patient. ESP32 microcontroller is connected to the GSM module as shown above. Initially, the ESP32 microcontroller is made to connect to the user's wifi. Whenever the GSM module is exposed to an open area, it will be connected to a number of satellites. Once it is connected, the user's location can be tracked by satellites and sent to the gsm module through transmitters and receivers. So this data will be available for ESP32 which can control it as per the requirements. Here we have set the boundary. So that whenever a patient goes beyond the boundary, it will notify the caregiver through the MY MQTT Application. This App is securely connected to the cloud where it receives all the data from ESP32.

GPS Tracking and Safe Zone Monitoring: An overview of the wristband's functionality, including its GPS tracking capabilities and the ability to monitor safe zones for the patient. This section explains how the device helps prevent wandering incidents and provides peace of mind to

caregivers.

Emergency Notifications: A description of how the wristband sends emergency notifications to caregivers if the patient ventures beyond the defined safe zones or requires immediate assistance. This section emphasizes the importance of timely response and enhanced patient safety.

Integration with Smart Spectacles: An explanation of the integration between the smart spectacles and the wristband, allowing seamless communication and data exchange. This section highlights the synergy between the two components to provide comprehensive support for Alzheimer's patients.




Figure 7.2: Block diagram of Wrist Band

7.6 Technological Considerations

Augmented Reality (AR) and its Role: An exploration of how augmented reality (AR) technology can be utilized in MemorEYES to provide visual cues, reminders, and additional information to Alzheimer's patients. This section discusses the potential benefits and challenges of implementing AR in the device.



Figure 7.3: The designed Wrist Band prototype

Virtual Reality (VR) for Cognitive Stimulation: A discussion on the integration of virtual reality (VR) technology in MemorEYES to offer cognitive stimulation exercises and immersive experiences for Alzheimer's patients. This section explores the potential of VR in improving cognitive functions and engagement.

GPS and Location-Based Services: An analysis of how GPS and location-based services play a crucial role in the device's functionality, including tracking the patient's location, setting safe zones, and providing real-time information to caregivers. This section highlights the importance of accurate and reliable location tracking.

Data Management and Security: A consideration of the data management and security measures implemented in MemorEYES to protect patient privacy and comply with regulatory requirements. This section discusses data encryption, secure storage, and user consent.

Usability and User Experience: A focus on the usability and user experience aspects of

MemorEYES, including user interface design, device comfort, intuitive operation, and accessibility features. This section emphasizes the importance of designing the device with the specific needs of Alzheimer's patients in mind.

7.7 Advantages of MemorEYES

MemorEYES offers several advantages for both Alzheimer's patients and their caretakers. These advantages include:

Users: The primary beneficiaries of MemorEYES are Alzheimer's patients and their caretakers. The device provides support, reminders, and safety features that enhance the patient's quality of life and reduce the caregiver's burden.

Portability and Customizability: MemorEYES is designed as a wearable device, ensuring portability and ease of use for patients. It can be customized to accommodate individual preferences, such as displaying specific photos or adjusting notification settings.

Risk Reduction for Caretakers: By assisting with medication reminders and GPS tracking, MemorEYES reduces the risk and stress faced by caretakers. They can have peace of mind knowing that the patient's medication schedule is being monitored, and they will be promptly alerted if the patient wanders beyond safe zones.

Long Durability: MemorEYES is built with durability in mind, considering the unique needs and challenges of Alzheimer's patients. The device is designed to withstand daily use and provide reliable performance over an extended period.

Enhanced Cognitive Function: The memory stimulation and emotional connection features of MemorEYES help improve cognitive function in early and mild-stage Alzheimer's patients. The display of pictures and reminders can trigger memories, enhance engagement, and promote mental well-being.

Improved Quality of Life: Overall, MemorEYES aims to improve the quality of life for Alzheimer's patients by providing support, connection, and safety. The device enables patients to maintain their independence and engagement with their surroundings, leading to a more fulfilling and empowered lifestyle.

7.8 Future Scope- Medical assisting kit

The MemorEYES device has significant potential for further advancements and expansion to enhance its functionality and address the needs of Alzheimer's patients. Some of the potential areas for future development include:

Navigation Capability: The smart glasses can be equipped with advanced navigation features to assist patients in finding their way back home or navigating unfamiliar environments. By integrating GPS technology and mapping services, MemorEYES can provide step-by-step directions, visual cues, and audio instructions to guide patients safely.

Mental Exercises and Engagement: New techniques can be incorporated into MemorEYES to engage patients in mental exercises that promote cognitive stimulation and memory retention. These exercises could include puzzles, memory games, and interactive activities specifically tailored to the needs of Alzheimer's patients.

Routine Reminders: MemorEYES can remind patients of important routines and activities such as drinking sufficient water, taking medication, and performing meditation or relaxation exercises. These reminders can be displayed on the glasses' interface or delivered through audio prompts, ensuring that patients adhere to healthy habits and self-care practices.

Health Monitoring Sensors: The inclusion of health monitoring sensors in MemorEYES can provide valuable insights into patients' well-being. Sensors such as oxygen level and heart rate monitors can help caregivers and healthcare professionals track vital signs and detect any abnormalities or health issues in real-time, enabling timely interventions.

Battery Life Improvements: Continued advancements in battery technology can extend the battery life of MemorEYES. Longer battery life would ensure that the device can be used throughout the day without frequent recharging, increasing convenience and usability for patients and caregivers.

7.9 Integration with Data analysis

MemorEYES, the smart wearable device designed for Alzheimer's patients, can greatly benefit from data analysis using Convolutional Neural Network (CNN) models. By analyzing the data collected from Alzheimer's patients, such as cognitive assessments, behavioral patterns, and physiological signals, MemorEYES can effectively detect early signs of the disease. This early detection is crucial as it allows for timely intervention and treatment, potentially slowing down the progression of Alzheimer's and improving patient outcomes. Through CNN analysis, MemorEYES can identify patterns and anomalies in the data that may indicate the presence of Alzheimer's disease, providing an early warning system for patients and caregivers.

One of the key advantages of using CNN models for data analysis in MemorEYES is the ability to personalize care planning. By analyzing the patient's data, such as memory tests, daily activities, and medication adherence, MemorEYES can create individual profiles that capture the unique characteristics and needs of each patient. This personalized approach enables healthcare professionals to develop tailored care plans that address specific challenges and support the cognitive and emotional well-being of Alzheimer's patients. By understanding the patient's patterns, preferences, and limitations, MemorEYES can provide targeted reminders, prompts, and assistance that cater to their specific requirements.

In addition to early detection and personalized care planning, data analysis using CNN models in MemorEYES offers the potential for continuous monitoring and remote patient assistance. By continuously analyzing the patient's data, including behavioral changes, sleep patterns, and physiological signals, MemorEYES can provide real-time insights into the patient's condition. This allows for ongoing treatment monitoring, ensuring that interventions and therapies are effective and can be adjusted as needed. Moreover, CNN-powered data analysis enables remote patient monitoring, where caregivers and healthcare professionals can access the patient's data and receive alerts or notifications regarding any concerning changes or risks. This remote monitoring capability provides convenience, peace of mind, and a more proactive approach to Alzheimer's care, reducing the need for frequent in-person visits and enabling timely interventions from a distance.

7.10 Conclusion- Medical assisting kits

Since ML algorithms provide data and analysis in real-time, doctors and other clinical experts can find problems much earlier. With machine learning in healthcare, many administrative tasks can be done automatically so that patients can get better care. In the early stages of researching new medicines, this technology can help. It is already ruled by research and development technologies like precision medicine and next-generation sequencing, which are used to find new ways to treat diseases that are hard to treat. By looking at a huge amount of patient data, doctors can find health problems before they turn into diseases. Clinical institutions can use ML to find strokes based on other health problems, check the health of the heart, and find other problems.

Using machine learning (ML) and deep learning (DL) models in clinical applications has a lot of potentials to change the way healthcare services are usually given. But different privacy and security problems need to be solved before these models can be used in clinical settings in a safe and reliable way. In this paper, we gave an overview of these problems by describing the ML pipeline in healthcare and pointing out the different places where it could be broken. We also talked about possible ways to make machine learning safe and private for security-critical applications like healthcare.

8 Future Scope-NeuroPredict

Expansion of the Dataset: To expand the dataset, it is recommended to incorporate additional sources of data such as PET scans and genetic data. PET scans provide valuable insights into brain activity and can contribute to a more comprehensive analysis of Alzheimer's disease. Genetic data, on the other hand, can help identify genetic risk factors and further enhance the accuracy of predictions. Collaborating with medical institutions will allow access to diverse datasets, ensuring a broader representation of the population and improving the model's generalizability. Additionally, the permission to access the ADNI dataset is a valuable opportunity as it contains a wealth of data related to Alzheimer's disease.

Long-term Monitoring and Predictive Analytics: Integrating wearables into the model's

data collection process can enable the tracking of physiological parameters over extended periods. This includes monitoring sleep patterns, heart rate variability, physical activity levels, and other relevant metrics. By incorporating these data points, the model can provide personalized interventions and detect early warning signs of cognitive decline or other related conditions. This long-term monitoring approach adds a dynamic aspect to the model and enables proactive interventions for better patient outcomes.

Accessibility and Scalability: To ensure accessibility and scalability, it is advisable to leverage cloud platforms. Cloud platforms provide the necessary infrastructure to process and store large volumes of data efficiently. They offer scalability options, allowing the model to handle increasing amounts of data without sacrificing performance. Additionally, investing in a powerful server infrastructure will support the computational demands of the model, ensuring fast and reliable processing. Implementing robust security measures is crucial to protect sensitive medical data and ensure data privacy compliance.

By expanding the dataset, incorporating textual data, integrating wearables, and ensuring accessibility and scalability, your machine learning model will be equipped with a more comprehensive and diverse set of information, leading to more accurate predictions and improved insights into Alzheimer's disease. These enhancements have the potential to advance the field of personalized medicine and contribute to early detection and intervention strategies.

9 Conclusions - NeuroPredict

Our project, NeuroPredict, successfully demonstrates the effectiveness of utilizing Convolutional Neural Networks (CNN) in predicting Alzheimer's disease based on MRI images. With its exceptional accuracy and efficiency, the CNN model employed in NeuroPredict provides a reliable tool for early detection and diagnosis of this devastating condition. By leveraging the power of CNN and combining it with Exploratory Data Analysis (EDA), our project offers a comprehensive solution for understanding the disease's underlying patterns and developing targeted interventions. Through the integration of EDA, we gain valuable insights into the demographic and clinical factors associated with Alzheimer's disease. This knowledge enables a deeper understanding of the

disease's progression, risk factors, and potential preventive measures. By identifying correlations and patterns within the data, we can develop proactive strategies that promote brain health, targeting vulnerable populations and providing tailored support to individuals already diagnosed with Alzheimer's disease.

The development of a hardware kit that incorporates the insights obtained from EDA further enhances our project's impact. This hardware kit not only aims to prevent Alzheimer's disease but also provides personalized support to individuals based on their specific needs. By analyzing various parameters such as gender, education, socioeconomic status, and age, we can offer targeted interventions that improve quality of life and potentially slow down the progression of the disease. Additionally, the user-friendly web application developed as part of NeuroPredict ensures accessibility and convenience for individuals seeking Alzheimer's prediction. By allowing users to directly upload their MRI images, we eliminate barriers and empower individuals to take control of their health. This web application serves as a vital tool in promoting early detection, encouraging timely medical interventions, and fostering informed decision-making. The significance of our project lies in its potential to make a substantial impact on the lives of individuals affected by Alzheimer's disease. By combining cutting-edge technology, data analysis, and personalized care, NeuroPredict paves the way for improved management, better patient outcomes, and a more proactive approach to Alzheimer's care. Furthermore, our project contributes to the broader research and understanding of neurodegenerative diseases, driving advancements in the field and shaping the future of healthcare. Furthermore, NeuroPredict offers a level of convenience that empowers individuals from all walks of life to upload their MRI images and receive accurate predictions. Recognizing the importance of accessibility and user-friendliness, our project has developed a user-friendly web application that simplifies the process for the common man.

Traditionally, obtaining a diagnosis for Alzheimer's disease has been a complex and time-consuming process, often involving multiple medical appointments and specialized facilities. However, with the NeuroPredict web application, individuals can conveniently upload their MRI images from the comfort of their own homes. This eliminates the need for physical visits to medical facilities and reduces the associated costs and time commitments. The web application has been designed with a focus on simplicity and ease of use, ensuring that individuals with little to no technical expertise can navigate the process effortlessly. The interface guides users through the steps required

to upload their MRI images, providing clear instructions and feedback along the way. By streamlining the process, we aim to remove any barriers that may prevent individuals from seeking early detection and timely medical interventions.

Once the MRI images are uploaded, the NeuroPredict CNN model swiftly analyzes the data and provides accurate predictions regarding the presence of Alzheimer's disease. The results are promptly delivered to the user, ensuring minimal waiting time and enabling individuals to take immediate action based on the outcome. This convenience and efficiency not only empower individuals but also foster a proactive approach to healthcare, where early detection and intervention can significantly impact the progression and management of the disease. By providing a convenient and accessible platform for individuals to upload their MRI images and obtain predictions, NeuroPredict revolutionizes the way Alzheimer's disease is diagnosed and managed. This novel approach empowers the common man to take control of their health and make informed decisions about their well-being. It eliminates geographical constraints and reduces the burden on healthcare systems, making Alzheimer's prediction more widely available to those who need it most.

In conclusion, NeuroPredict showcases the power of CNN models in predicting Alzheimer's disease, the insights gained from EDA, and the development of innovative hardware and web applications. By focusing on early detection, proactive prevention, and personalized care, our project addresses the challenges posed by Alzheimer's disease and offers a comprehensive solution. As we continue to refine and expand our project, we strive to make a lasting impact on individuals affected by Alzheimer's disease and contribute to the global fight against neurodegenerative disorders.

REFERENCES

- [1]Zellefrow, B., Shanaei, S., Salam, R., Nasan, A. E., Elzabadani, H. (2017). Halcyon-Assistive Technology for Alzheimer'S Patients. In Recent Advances in Information Systems and Technologies: Volume 2 5 (pp. 375-383). Springer International Publishing.
- [2]Lehrer, P., Kaur, K., Sharma, A., Shah, K., Huseby, R., Bhavsar, J., Zhang, Y. (2020). Heart rate variability biofeedback improves emotional and physical health and performance: A systematic review and meta analysis. *Applied psychophysiology and biofeedback*, 45, 109-129.
- [3]Salehi, W., Gupta, G., Bhatia, S., Koundal, D., Mashat, A., Belay, A. (2022). IoT-Based Wearable Devices for Patients Suffering from Alzheimer Disease. *Contrast Media Molecular Imaging*, 2022.
- [4]Gacem, M. A., Alghlayini, S., Shehieb, W., Saeed, M., Ghazal, A., Mir, M. (2019, December). Smart Assistive Glasses for Alzheimer's Patients. In 2019 IEEE International Symposium on Signal Processing and Information Technology (ISSPIT) (pp. 1-5). IEEE.
- [5]O'Bryant, S. E., Edwards, M., Johnson, L., Hall, J., Villarreal, A. E., Britton, G. B.,... Graft-Radford, N. R. (2016). Ablood screening test for Alzheimer's disease. *Alzheimer's Dementia: Diagnosis, Assessment Disease Monitoring*, 3, 8 3 - 90. doi:10.1016/j.dadm.2016.06.004
- [6]Ebrahimi-Ghahnavieh, A., Luo, S., Chiong, R. (2019). Transfer Learning for Alzheimer's Disease Detection on MRI Images. 2019 IEEE International Conference on Industry 4.0, Artificial Intelligence, and Communications Technology (IAICT). do:10.1109/iciaict.2019.8784845
- [7] Abdulkadir, A., Mortamet, B., Vemuri, P., Jack Jr., C.R., Krueger, G., Kl€oppel, S., Alzheimer's Disease Neuroimaging Initiative, 2011. Effects of hardware heterogeneity on the performance of SVM Alzheimer's disease classifier. *NeuroImage* 58 (3), 785-792.
- [8] Kourtis, .L C., Regele, O. B., Wright, J. M., Jones, G. B. (2019). Digital biomarkers for Alzheimer's disease: the mobile/wearable devices opportunity. *Npj Digital Medicine*, 2(1). do:10.1038/41746-019- 0084-2
- [9] Agüera-Ortiz, A., et al. "A Survey of Alzheimer's Disease Early Diagnosis Methods for Cognitive Assessment." *Frontiers in Aging Neuroscience* 12, no. 185 (2020): 1-13. doi:10.3389/fnagi.2020.00185.
- Figure 1.